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Instructions for Use

MD

Intended Purpose
This non-sterile, manual, hand-held device is intended for dispensing two-component dental materials (e.g. elastomeric impression materials, dental cements) from a cartridge directly into the patient’s mouth or into another device (e.g. impression tray). This is a reusable device.

Indications for use
Indicated for use with dental material such as elastomeric impression materials and dental cements for cementing, resin-based composite materials for core build-ups or Bis-acryl materials for temporary restorations.


Performance features
The performance features of the product meet the requirements of the intended use.

Contraindications
There are no known contraindications.


Patient target group
Persons who are treated during a dental procedure.

Intended users
This medical device should only be used by a professionally trained dental practitioner.

Application
1. Push the small lever (1) in the back of the AUTOMIX Dispenser and retract the slide (2) completely. Lift the lever (3) on the top of the Dispenser and insert the cartridge.

2. Remove the cartridge cap (4) by turning it 90° counter-clockwise. Attach a mixing cannula. Make sure that the guidance of the AUTOMIX-cartridge is aligned with that of the mixing cannula and turn the cannula 90° clockwise until it locks in position.

3. The material is dispensed through the mixing cannula by pulling the trigger (5) of the AUTOMIX Dispenser and can be applied directly.

Removal of the cartridge:
Push up lever (1) and retract the slide. Lift lever (3) on the top of the Dispenser and remove the cartridge.

Store used cartridge with fixed used mixing cannula. It serves as a cap.

General Information
This dispenser is delivered non-sterile.

Direction for processing according EN ISO 17664

Reprocessing
Disassembly

Limitations on processing

- Any dispenser showing visible damage or wear must be exchanged.
- End of life of a reusable dispenser is normally determined by wear and damage due to use.
- Use-life of the dispenser is exceeded if visual inspection shows unacceptable deterioration like wear, discoloration, crackles, etc.

Initial treatment at the point of use

- To prevent any contamination from drying it is recommended that reprocessing takes place not later than 2 hours after use.
- Keep the dispenser moist after use to prevent soil from drying on them.
- Contaminated dispenser should be transported to the area for cleaning in a clean and closed box/container that avoids spreading of the contamination.

Manual cleaning

- Use a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) to clean all surfaces of the plunger for 30 seconds. Replace wipe, as necessary, if it becomes soiled.

- In order to remove residue on the plunger, use a soft bristle brush or a toothbrush and brush for 30 seconds. If necessary, increase time until no visible soil is observed.
- Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces of the plunger two (2) times to remove residual detergent agent. Use additional wipes as needed.
- Place a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) through the gaps of the plunger and move the wipe back and forth for 30 seconds in the gaps to clean the plunger. Replace wipe, as necessary, if it becomes soiled.
- In order to remove residue in the gaps of the plunger, use a soft bristle brush or a toothbrush and brush for 30 seconds. If necessary, increase time until no visible soil is observed.
- Use a clean, lint-free wipe, dampened with purified water to wipe all gaps of the plunger two (2) times to remove residual detergent agent. Use additional clean, dampened wipes as needed.
- Use a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) to clean all surfaces of the dispenser for 30 seconds. Replace wipe, as necessary, if it becomes soiled.
- In order to remove residue in the grooves of the dispenser, use a soft bristle brush or a toothbrush and brush for 140 seconds. If necessary, increase time until no visible soil is observed.
- Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces of the dispenser two (2) times to remove residual detergent agent. Use additional wipes as needed.
- Place a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) through the gaps of the dispenser and move the wipe back and forth for 30 seconds to clean the gaps. Replace wipe, as necessary, if it becomes soiled.
- In order to remove residue in the gaps, use a soft bristle brush or a toothbrush and brush for 30 seconds. If necessary, increase time until no visible soil is observed.
- Use a clean, lint-free wipe, dampened with purified water to wipe all gaps of the dispenser two (2) times to remove residual detergent agent. Use additional clean, dampened wipes as needed.

Manual disinfection

- Use a clean, lint-free wipe, dampened with a disinfectant containing isopropyl alcohol (e.g. CaviWipes 17%) to allow all surfaces and gaps of the plunger and dispenser to remain visibly wet for a minimum of five (5) minutes. Use additional wipes to ensure all surfaces and gaps remain wet for the five (5) minute contact time.
- Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces and gaps of the plunger and dispenser for at least two (2) times to remove residual disinfectant. Use additional wipes as needed.

Drying

- Allow the plunger and dispenser to dry until visibly dry.

Maintenance, Inspection and Testing
Inspect the dispenser and plunger under normal lighting:

- Any dispenser or plunger showing visible damage or wear must be exchanged.
- Repeat manual cleaning and inspection if the dispenser or plunger are not visibly clean.

Packaging
Reassemble the dispenser and plunger.

Sterilization
Dispenser are non-sterile medical device. Sterilization is not necessary.

Disposal
Disposal of the product according to local authority regulations.

Reporting obligation
Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Warranty
P. L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P. L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and P. L. Superior Dental Materials GmbH’s sole obligation shall be repair or replacement of the P. L. Superior Dental Materials GmbH product.

Limitation of Liability
Except where prohibited by law, P. L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

Lang Dental Mfg. Co., Inc.
175 Messner Dr
Wheeling, IL 60090
USA