

Splintline[®] Acrylic

Intended Purpose

Fabrication, by a dental professional, of fast-curing provisional crowns and bridges intended to remain in the mouth for less than 30 days.

Other Materials / Equipment Needed

1. Silicone lubricant
2. Aquapres[®] (pressure curing unit), optional
3. Acrylic bur or stone
4. Jet[™] Adjusters (acrylic stains), optional
5. Jet[™] Seal (acrylic sealant), optional

Instructions for Use



For consistent results, always shake liquid before use, as catalyst may have separated from monomer.

1. In the impression, fill the tooth that was extracted or prepared with Splintline[®] as follows: Moisten surface with Splintline[®] liquid. Spray on Splintline[®] powder. Alternate liquid and powder until impression is filled. For best results, begin and end with liquid.
2. While allowing acrylic to reach a paste like consistency, coat the prepared tooth and surrounding tissue with a silicone lubricant to allow easy removal of acrylic from prep and to protect the teeth and surrounding tissue.
3. Seat impression in mouth and hold in place for three minutes until material has reached a rubber like state. Remove and trim excess acrylic. Do not over trim.
4. Tease the Splintline[®] provisional from the impression and replace into patient's mouth. Have patient bite into occlusion for 60 seconds. CAUTION! DO NOT allow the material to get hot in the mouth!
5. Remove provisional from patient's mouth. Allow an additional 10 minutes for curing, either on bench or, for better results in Aquapres[®] at 30 psi.

6. Check fit. Trim excess material with acrylic bur or stone. Grind to a smooth finish.
7. Modify shade with Jet™ Adjusters if necessary. Polish using Jet™ Seal.
8. Clean patient's teeth of any residual silicone lube. Place provisional into mouth and secure with suitable temporary cement.

Handling Precautions

Contains ethyl methacrylate

Warning. Flammable liquid and vapor. Keep away from heat/sparks/open flames/hot surface. – No smoking.

Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation.

Wear protective gloves/protective clothing/eye protection/face protection.

Keep material in original container tightly closed in upright position. Store in a well-ventilated place. Keep cool.

Disposal

Dispose of contents and container in accordance with local, state, and national legislation.

Warranty

Manufacturer warrants this product to be free from manufacturing defects for one year from date of purchase only from an authorized reseller. Warranty is void if product is not used for intended purpose, not used as instructed or not stored in original container.

Lang Dental Mfg. Co., Inc.
175 Messner Dr.
Wheeling, IL 60090 USA
www.langdental.com
Ph. +1 847-215-6622
Fax +1 847-215-6678
Made in USA

LANG[®]



MediMark® Europe
11 rue E. Zola
38100 Grenoble
France



Note: Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

©2021 Lang Dental Mfg. Co., Inc. LANG, Splintline and AQUAPRES are the registered trademarks of Lang Dental Mfg. Co., Inc. JET is the trademark of Lang Dental Mfg. Co., Inc. All rights reserved.

LDD 20 R. 7 2021-02

