

Denture Duplicator Flask

Intended Purpose

Lab Processing of:

- Duplicate Denture from a patient's original denture
- Surgical Guide/Implant Stent

DENTURE DUPLICATION

Materials / Equipment Needed

Silicone Spray

Alginate

Spatula

Jet™ Liquid

Dropper

Jet™ Tooth Shade Powder

Bur, pumice, lab polish

Jet™ Denture Repair Powder

Measuring cup


Aquapres®/pressure curing unit
(optional)

Jet™ Adjusters (optional)

Jet™ Seal (optional)

Instructions for Use

ALWAYS SHAKE JET™ LIQUID BEFORE USE.



- Spray inside of flask with silicone spray to facilitate post procedure clean-up.
- Fill half of flask with alginate.



- Place existing denture into alginate. Teeth should be perpendicular to the bottom of the flask.
- With wet fingers, smooth alginate from borders of denture to rim of flask, removing excess alginate.



- After alginate is set, trim excess that covers the denture.



- Spray denture and surrounding alginate with silicone spray.
- Fill other half of flask and palate side of denture with alginate.
- With wet fingers, smooth alginate to ensure complete coverage, working quickly before alginate sets.



- Close flask with steady pressure. Secure thumb screw and allow alginate to set.
- Trim excess alginate.
- Open flask and remove denture. If there are any gross defects in the impression, make corrections now or retake impression.



- Select appropriate Jet™ Tooth Shade powder. Fill teeth in impression as follows: Moisten surface with Jet™ liquid. Spray on Jet™ powder. Alternate liquid and powder until teeth are filled to gingival border. When acrylic appears dull, use small spatula to push acrylic down to redefine gingival borders of all teeth.



- Cure for 15 minutes, either on bench or, for best results, in Aquapres® at 30 psi. (Curing under pressure eliminates porosity, making a strong, dense, aesthetically superior denture.)



- With bur, trim teeth interproximally. Place trimmed teeth back into impression.
- Place 30 cc of Jet™ Denture Repair powder into measuring cup. Add 15 cc of Jet™ liquid and mix for 10 seconds.
- Pour 75% of mixture into lower half of flask. Pour remaining acrylic over palate on upper half of flask.
- With slow and steady pressure, close and secure flask. (Excess speed or pressure may cause distortion.)
- Remove excess acrylic flashing from outside of flask, as hardened acrylic will inhibit movement of hinge and thumb screw.
- Cure for 15 minutes as previously instructed.



- Remove denture from flask, clearing away all excess acrylic.
- Grind to a smooth finish.
- Polish with pumice, then high shine.



- Use Jet™ Adjusters for polishing and characterizing teeth.
- Finish with Jet™ Seal to fill in any surface porosity and leave a shiny, finished surface.

CLEAR SURGICAL GUIDE / IMPLANT STENT

Materials / Equipment Needed

Silicone Spray

Alginate

Spatula

Jet™ Liquid

Dropper

Measuring cup

Ortho-Jet™ Powder
(Clear Guide)

Jet™ XR Opaque

(Radiopaque Appliance)

Jet™ XR Shadow

(Radiopaque Appliance)

Bur, pumice, lab polish

Aquapres®/pressure curing unit
(optional)

Instructions for Use – Clear Guide

Follow Steps 1-5 as shown for Denture Duplication.

6. Place 30 cc of Ortho-Jet™ powder into a measuring cup. Add 15 cc of Jet™ liquid and mix for 10 seconds. Pour 75% of mixture into the impression of the denture in the lower half of flask. Ensure that acrylic flows completely into all incisal and occlusal spaces without trapping any air.

Pour the rest of the mixture over the palate on the upper half of the flask.

7. With slow and steady pressure close and secure flask. (Excess speed or pressure may cause distortion.) Cure for 20 minutes either on bench or, for best results, in Aquapres® at 30 psi.

Follow Step 9 for the Denture Duplication to complete procedure.

Instructions For Use — Radiopaque Appliance

Follow steps 1-5 as shown for Denture Duplication.

6. Follow instructions for Denture Duplication, replacing the Jet™ Tooth Shade powder with the Jet™ XR Opaque powder.
7. Same as for Denture Duplication.
8. Follow instructions for Denture Duplication, replacing the Jet™ Denture Repair powder with the Jet™ XR Shadow powder.

Follow Step 9 for Denture Duplication to complete procedure.

Warranty

Manufacturer warrants this product to be free from manufacturing defects for one year from date of purchase. Warranty is void if product is not used for intended purpose or not used as instructed.

Lang Dental Mfg. Co., Inc.
175 Messner Drive
Wheeling, IL 60090 USA
www.LangDental.com

LANG®

Ph +1 847-215-6622
Fax +1 847-215-6678



©2020 Lang Dental Mfg. Co., Inc. LANG and AQUAPRES are the registered trademarks of Lang Dental Mfg. Co., Inc. JET, JET XR and ORTHO-JET are the trademarks of Lang Dental Mfg. Co., Inc. All rights reserved.

LDD 0398 R. 9 2020-01