



Sophisticated, innovative, and unique

Aquapres™ by Lang Dental — While visiting the exhibits at this year's



Aquapres™ by Lang Dental

Hinman Dental Meeting, I was talking to Dave Lang, the owner of Lang Dental Company. He asked me if I had a Lang Aquapres,™ and I said that I thought every dentist had one. He said, “No they don’t.” I could not believe it because I have used one for years. I asked Dave if he had raised the price. He said he had not, and that the Aquapres™ is still very inexpensive. I said maybe the problem is that it never wears out. Just in case you are one of those who do not know

about the Lang Aquapres™, it will make anything you do with acrylic resins — such as provisionals or removables — better. One of the easiest and most trouble-free methods to obtaining stronger, denser, and porosity-free acrylic work is to use the Lang Aquapres™ hydraulic pressure curing unit. It works on relines, repairing acrylic on partials. Aquapres™ is safe to use because it requires no compressed air or electricity. It operates safely under hydraulic pressure, and is easily applied with a twist of the wrist. The complete fabrication system consistently generates dependable, high-quality results by eliminating porosity in methacrylate acrylics while achieving long-lasting color. Equally appealing is the low cost of this unit. Small and compact, Aquapres™ can also easily be stored in your lab. It is portable, not only within your office, but also to remote clinics such as nursing homes. High quality, dependable, and color-stable, acrylic work is easily, conveniently, and quickly accomplished in your office, your office lab, or chairside. Aquapres™ is a problem solver that every practice will find helpful. Don't be left out in the cold without the Lang Aquapres™. Order one via any dental dealer, or visit www.langdentalmanufacturing.com for more information. **DE**

