

Effective Date 08-Jun -2023

Version 4

<b>1. IDENTIFICATION</b>
--------------------------

**Product Identifier****Product Name** RESINCAP POWDER**Other means of identification****SDS#** 003**Product Code** 2310, 2312, 2350, 2370, 2380, 2390**Recommended use of the chemical and restrictions on use****Recommended Use** Fabrication of patterns and copings**Details of the supplier of the safety data sheet****Supplier Address** Lang Dental Mfg. Co., Inc.  
175 Messner Dr.  
Wheeling, IL 60090  
USA**Emergency telephone number****Company Phone Number** +1 847-215-6622**Emergency Telephone (INFOTRAC)** +1 352-323-3500 (International)  
800-535-5053 (North America)**Authorized European Representative**MediMark® Europe SARL  
11 rue Emile Zola. BP 2332  
38033 Grenoble Cedex 2  
France  
Tel: +33 (4) 76 86 43 22  
Fax: +33 (4) 76 17 19 82  
Email: info@medimark-europe.com

<b>2. HAZARDS IDENTIFICATION</b>
----------------------------------

**Classification of the substance or mixture**

Eye Damage/Irritation	2B
Skin sensitization	1

**Signal word** Warning**Hazard statements** H317 May cause an allergic skin reaction.  
H320 Causes eye irritation.**Physical State** Powder      **Appearance** Fine, red      **Odor** Faint odor in bulk**Precautionary Statements – Prevention**

- P202 Do not handle until all safety precautions have been read and understood.
- P261 Avoid breathing dust, fumes, gas, mist, vapors or spray.
- P264 Wash hands and exposed skin thoroughly after handling.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.

### Precautionary Statements – Response

P363 Wash contaminated clothing before reuse.  
P302+P352 IF ON SKIN: Wash with soap and water.  
P305+P351 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy  
+P338 to do – continue rinsing.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
P337+P313 Get medical advice/attention.

### Precautionary Statements – Disposal

P501 Dispose of contents/container to an approved waste disposal plant.

## 3. COMPOSITION / INFORMATION ON INGREDIENTS

### Classification of the substance or mixture

Chemical Name	CAS No	Weight - %	Trade Secret
Polymer	9011-14-7	< 80	*

\*Specific chemical weight has been withheld as a trade secret.

## 4. FIRST AID MEASURES

### First aid measures

<b>Inhalation</b>	Remove victim to fresh air and keep at rest in a position comfortable for breathing. Get medical attention if discomfort persists.
<b>Eye contact</b>	Rinse immediately with plenty of water, including under the eyelids, for at least 15 minutes. If irritation persists, call a physician immediately.
<b>Ingestion</b>	Do NOT induce vomiting. Drink plenty of water or milk immediately. If vomiting, continue to offer water or milk. Never give anything by mouth to an unconscious person. Call a physician or poison control center immediately and provide an estimate of when and how much material was ingested.
<b>Skin Contact</b>	Wash with soap and water. If irritation, redness or swelling persists, call a physician immediately. Take off contaminated clothing and wash before reuse.

### Most important symptoms and effects, both acute and delayed

**Symptoms** No information available.

### Indication of any immediate medical attention and special treatment needed

**Note to physicians** Treat symptomatically.

## 5. FIRE-FIGHTING MEASURES

### Extinguishing Media

**Suitable:** Chemical (alcohol-resistant) foam, dry chemical or carbon dioxide.

**Unsuitable:** Water may not be effective in extinguishing this fire.

### Specific hazards arising from the chemical

For bulk size: Polymers are combustible dusts, care should be taken to avoid creating explosive concentrations in the air. Follow grounding and bonding procedures.

### Special Fire Fighting Procedures:

Avoid extinguishing methods, which may generate dust clouds. Water stream can disperse dust into air producing a fire hazard and possible explosion hazard if exposed to ignition source. Firefighters should wear self-contained breathing apparatus.

### Protective equipment and precautions for firefighters

For bulk size: Polymer dust is combustible. The explosive limits of the polymer particles suspended in air are approximately those of coal dust. Polymers are sensitive to static charge, follow grounding and bonding procedures. Polymers are not sensitive to mechanical impacts.

## 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

**Personal precautions** Use personal protective equipment as required. Remove any contaminated clothing and wash thoroughly before reuse.

### Methods and material for containment and clean-up

**Method for containment** Prevent further leakage or spillage if safe to do so. May contaminate water supply.

**Method for clean-up** Sweep up to avoid slipping hazard. Keep airborne particulates at a minimum when cleaning up spills. Clean up in accordance with all applicable regulations. Wash all affected areas with plenty of warm water and soap.

## 7. HANDLING AND STORAGE

### Precautions for safe handling

**Advice on safe handling** Use only in well-ventilated areas. Avoid contact with skin, eyes or clothing. Avoid breathing dust or fume. Handle in accordance with good industrial hygiene and safety practice. Wash thoroughly after handling. Do not eat, drink or smoke when using this product.

### Conditions for safe storage, including any incompatibilities

**Storage Conditions** Keep container tightly closed to prevent water absorption and contamination. Store in a dry, cool and well-ventilated place away from direct sunlight or other sources of light or intense heat. The temperature should remain at normal temperature not to exceed 30°C. Storing above recommended temperature will cause product performance issues. Store in accordance with federal, state and local regulations. Observe all label precautions until the container is cleaned, reconditioned, or destroyed.

**Packaging materials** Keep in original container.

**Incompatible materials** Strong oxidizers, strong oxidizing agents

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Appropriate engineering controls** Apply technical measures to comply with the occupational exposure limits. When working with large quantities of product, provide adequate ventilation (e.g. local exhaust ventilation, fans). Ensure that an eyewash station, sink or washbasin is available in case of exposure to eyes. Use good local exhaust at processing equipment, including buffers, sanders, grinders and polishers.

### Individual protection measures, such as personal protective equipment

**Eye / face protection** Depending on the use of this product, safety glasses or goggles may be worn. If necessary, refer to US OSHA 29 CFR SS1910.133, Canadian standards or the European Standard EN 166. Ensure that an eyewash station, sink or washbasin is available in case of exposure to eyes.

**Skin and body protection** If anticipated that prolonged and repeated skin contact will occur during use of this product, wear gloves for routine industrial use. If necessary, refer to US OSHA 29 CFR SS1910.138 or the appropriate standards of Canada or the EC member states.

**Respiratory protection** No special respiratory protection is required under typical circumstances of use or handling. If necessary, use only respiratory protection authorized per US OSHA requirement in 29 CFR SS 1910.134, or applicable US state regulations, or the appropriate standards of Canada, its provinces, or the EC member states. VENTILATION: Local exhaust at processing equipment.

**General hygiene considerations** Handle in accordance with good industrial hygiene and safety practice. Wash thoroughly after handling. Food, beverages and tobacco products should not be carried, stored, or consumed where this material is in use. Wash hands thoroughly before eating, drinking, or smoking.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

<b>Physical state</b>	Powder	<b>Odor</b>	Faint odor in bulk
<b>Appearance</b>	Fine		
<b>Color</b>	Red		

<b>Property</b>	<b>Values</b>
<b>Boiling point / boiling range</b>	Not applicable
<b>Auto ignition temperature</b>	N/A
<b>Flash point</b>	N/A
<b>Flammability limits in air</b>	
<b>Upper flammability limit</b>	Not applicable
<b>Lower flammability limit</b>	Not applicable
<b>Specific gravity</b>	40000.012

## 10. STABILITY AND REACTIVITY

<b>Reactivity</b>	Not reactive under normal conditions
<b>Chemical stability</b>	Stable under recommended storage conditions
<b>Possibility of hazardous reactions</b>	None under normal processing
<b>Hazardous polymerization</b>	Does not occur.
<b>Incompatible materials</b>	Strong oxidizing agents
<b>Hazardous decomposition products</b>	Methacrylate monomer, oxides of carbon when burned

## 11. TOXICOLOGICAL INFORMATION

### Information on likely routes of exposures

<b>Product information</b>	This product has not been tested on animals to obtain toxicology data.
<b>Inhalation</b>	Not expected to be an inhalation hazard under normal conditions of intended use
<b>Eye contact</b>	Avoid contact with eyes.
<b>Skin contact</b>	Avoid contact with skin.
<b>Ingestion</b>	Do not swallow.

### Information on physical, chemical and toxicological effects

**Symptoms** No information available.

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

No information available

## 12. ECOLOGICAL INFORMATION

<b>Component Ecotoxicity</b>	There is no specific data available for this product; however, very large releases may be harmful or fatal to overexposed aquatic life.
------------------------------	---

## 13. DISPOSAL CONSIDERATIONS

### Waste treatment methods

<b>Disposal of wastes</b>	Follow all local and national government regulations in disposing material or contaminated packaging.
<b>Contaminated Packaging</b>	For bulk only: Reuse of empty drums or containers is not recommended. Employees should be advised of the potential hazards due to residual material associated with empty containers. Dispose of all empty containers properly in accordance with federal, state and local regulations.

## 14. TRANSPORTATION INFORMATION

<u>DOT</u>	Not regulated
<u>IATA</u>	Not regulated
<u>IMDG</u>	Not regulated

## 15. REGULATORY INFORMATION

### International Inventories

EINECS This product is a medical device and not subject to chemical notification requirements.

### US Federal Regulations

TSCA This product is an FDA regulated device and not subject to TSCA regulations.

SARA 302 – Extremely Hazardous – Not Listed

SARA 311 / 312 Hazard Categories – This product is an FDA regulated device and not subject to reporting requirements.

SARA 313 – Toxic Chemicals – Not established.

### US State Regulations

California Proposition 65 – Warning. This product contains the following chemicals which are listed by the State of California as carcinogenic or a reproductive toxin: None

US State Right-to-know Regulations Not established.

## 16. OTHER INFORMATION

HMIS	Health Hazards	Flammability	Physical Hazards
	0	0	0

**Effective Date** 08-Jun-2023

### Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release. It is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet