

LANG

The Industry Leader
in Acrylics

Tooth Shade Conversions

Lang Jet Tooth Shades are available in the Nue Hue* and Vita** shades. Use the chart below for your reference.

Nue Hue*	Lang	Vita**
59	59	B1
60	60	A1
61	61	B2
62	62	A2
65	65	D3
66	66	A3
67	67	B3
69	69	C2
77	77	B4
81	81	A3.5
87	87	A4
	C3	C3
	C4	C4
	D4	D4
	Light Incisal	
	Dark Incisal	
	Clear	
	OM1	OM1
	OM2	OM2

* Registered trademark of Dentsply, Int.

** Registered trademark of Vita Zahnfabric.

Shades refer to body shade only. When looking at the tooth, ignore the incisal 1/3 of the tooth.

Lang Dental Mfg. Co., Inc.

P.O. Box 969, 175 Messner Drive., Wheeling, Illinois 60090 USA
847-215-6622 • 800-222-5264 • Fax: 847-215-6678 • www.langdental.com